



908-707-8680

topsafety.com

AED Wall Cabinet: Surface Mount with Alarm, Security Enabled



This well-constructed AED wall cabinet is designed to store your AED in an easily accessible location. The AED sits inside the cabinet with a transparent door for better visibility. The AED can be stored inside the cabinet with or without a carry bag.

Sizing: Internal (H/W/D): 14.25 in (36.20 cm) X 14.25 in (36.20 cm) X 7 in (17.78 cm)

External (H/W/D): 17.5 in (44.45 cm) X 17.5 in (44.45 cm) X 7 in (17.78 cm)

Protrusion from the wall: 7 in (17.78 cm)

Material: Stainless steel

Security: Audible alarm triggered when the door is opened. The alarm can also be wired into a building's security system and can be enabled or disabled by provided key.

Included: All hardware and two key





